

ITN Flow Cytometry SOP

Commonly used panels (address questions to biologists@immunetolerance.org):

Cross-trial T cell panel

Target	Clone
CD3	OKT3
CD4	RPA-T4
CD8	SK1
CCR7	G043H7
CD45RA	HI100
CD45RO	UCHL1
CD127	A019D5
FoxP3	259D/C7
CXCR5	
PD-1	EH12.2H7
CD57	
KLRG-1	
Eomes	
T-bet	
CD56	NCAM16.2
TIGIT	MBSA43

B Cell Panel

Target	Clone
CD20	2H7
CD21	B-ly4
CD24	ML5
CD38	HB-7
CD27	O323
CD19	HIB19
IgM	MHM-88
IgD	IA6-2
CD3 (dump)	UCHT1
CD138	
CD45	
CD34	
CD14 (dump)	